



Specialized Equipment for Specialty Healthcare



MTI's procedure chairs and specialized exam room medical equipment offer a higher level of comfort and luxury with better caregiver access to properly examine and treat patients. Our ADA-compliant chairs feature tapered and thin backs giving better medical staff access to patients while also offering lower seat entry heights for easier access for disabled patients. These chairs come with many options including powered swivel bases, floating patient arms, hand/forearm procedure arms, and convenient headrest options. MTI's intuitive chairs are technology-rich and built to fit the needs of the specialist. Find out today how your practice can lower its operating costs, increase practice safety, and improve patient outcomes with MTI equipment.

MTI Chairs/Tables



	463/464 Exam/Procedure Chairs	829/830 Procedure Chairs	840 Stretcher Bed	528/529 Podiatry/Wound Care Chairs	529W Bariatric Podiatry/Wound Care Chair	550 Podiatry/Wound Care Chair	721 Surgery Chair
Seat Height Range	18.75" - 37"	18.75" - 39"	18.75" - 37"	18.75" - 39" (529), 21" fixed (528)	18.75" - 39"	18.75" - 39"	19" - 35"
Seat Width	21"	22" (829) or 25" (830)	25"	24"	30"	25"	21.5"
Patient Capacity	725 lbs	725 lbs	725 lbs	725 lbs	850 lbs	725 lbs	650 lbs
Programmable Positions	4	4	4	4	4	4	3
Base Options	Swivel	Standard, swivel or mobile	Mobile steering only	Standard or swivel	Standard	Standard, swivel or mobile	Standard, swivel or mobile
Controls	EasyView™ chair back mounted membrane switches	Includes both hand and foot controls	Includes sealed membrane hand control	4-function (529), 3-function (528), programmable foot control, optional hand control	4-function programmable foot control, optional hand control	4-function programmable foot control, optional hand control	Choice of either hand or foot control
Legrest Design	Leg rest rises as back is reclined	Independent footrest & leg extension	Independent footrest & leg extension	Fixed angle w/11" motorized extension	Fixed angle w/11" motorized extension	Independent footrest & leg extension	Comfortable fixed angle
Standard Armrests	Removable slide-back patient arms	Multiple optional arm configurations	Multiple optional arm configurations	Removable slide-back patient arms	Removable slide-back patient arms	Removable slide-back patient arms	Multiple optional arm configurations


Why MTI Chairs/Tables are Your Best Option



Ultra quiet DC motors are 30% faster and at 725 lb rating, our latest chairs have an unmatched patient rating



The soft and sleek design defines luxury and provides superior patient access with multiple spline locations to attach armrests and accessories..



Articulating headrests adjust and lock firmly into ideal positions in just seconds.



Supported by a welded steel frame with maintenance-free joints and unsurpassed strength for better sturdiness under greater loads.



A fundamental advantage to MTI's products is the tapered back section which allows easier access to the patient's head and neck area with less strain on the caregiver.



ADA Compliant chairs with an 18.75" entry height feature automatic power footrest extensions with extend/retract buttons for independent movements.

ADA Compliance

ADA Regulations Protect Patients and Care Providers by Avoiding Injuries and ensure accessibility to, and usability of diagnostic equipment by patients with disabilities.

MTI makes the only ADA-compliant M302 chairs on the market today which can help reduce patient and provider injury and distress with lower entry heights, advanced control features and functionalities. Elevate your patient and caregiver experiences with increased safety and comfort with MTI's technology-rich, intuitive, ADA compliant chairs.



Other MTI Products



Cabinets



Instrument Tables



Stools



Side Chairs



Lights

Types of Cabinets: 9
Countertop Height: 33.5 - 41.5"
Width: 25-35.5"
Drawers: 2-6
Finish: laminate or stainless steel
Counter: laminate or solid surface

Locking casters
Counter Finish: laminate or solid surface
Counter Size: 23.75" x 39" or 24.75" x 40"
17.5" vertical adjustment
Optional pull-out tray

Types of Stools: 12
Controls: hand or foot
Height: 18.25-28.75"
Seat Diameter: 16.5"
Backrest Options: none, ergo, curved body or twin arms
Weight Capacity: 225-350 lbs.

Types of Chairs: 2
Seat Height: 17.75"
Seat Width: 18.4"
Backrest: ergonomic
Weight Capacity: 425 lbs.
Arms: with or without
Adjustable leveling feet

Types of Lights: 3
Model: Clarus™ or Lumerus™ LED
Lux: 21,500 - 33,400
Life: 40,000 - 70,000 hrs
CRI: ≥ 75
Style: caster base, floor stand, or chair-mounted